



Recombinant *Macaca fascicularis* Inorganic pyrophosphatase (PPA1)

Product Code	CSB-EP680994MOV-B
Abbreviation	PPA1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4R543
Product Type	Recombinant Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	>85% (SDS-PAGE)
Sequence	SGFSTEERA APFSLEYRVF LKNEKGQYIS PFHDIPIYAD KDVFHMVVEV PRWSNAKMEI ATKDPLNPIK QDVKKGKRLRY VANLFPYKGY IWNYGAIPT WEDPGHNDKH TGCCGDNDPI DVCEIGSKVC ARGEIIGVKV LGILAMIDEG ETDWKVIAIN VDDPDAANYN DINDVKRLKP GYLEATVDWF RRYKVPDGKP ENEFAFNAEF KDKDFAIDII KSTHDHWKAL VTKKTNGKGI SCMNTTVSES PFKCDPDAAR AIVDALPPPC ESACTVPTDV DKWFFHHQKN
Source	<i>E.coli</i>
Target Names	PPA1
Protein Names	Recommended name: Inorganic pyrophosphatase EC= 3.6.1.1 Alternative name(s): Pyrophosphate phospho-hydrolase Short name= PPase
Expression Region	2-289
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.