



Recombinant Drosophila melanogaster Succinate dehydrogenase assembly factor 2-A, mitochondrial (CG14757)

Product Code	CSB-EP681172DLU-B
Abbreviation	CG14757
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4V5I9
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	SNMKDSP PPPPPLASTF NDVIVDYEDP DYLPPLPEYYPV RPNEPLETRK QRLLYQSRKR GMLENDLLLS TFAAKHLQNF SAEQTAQYDQ LINGVSNDWD IYYWATEVKP TPKEYDTEIM GLLKEHVKNA ERVTRLRQPD LNI
Source	E.coli
Target Names	CG14757
Protein Names	Recommended name: Succinate dehydrogenase assembly factor 2-A, mitochondrial Short name= SDH assembly factor 2-A Alternative name(s): Succinate dehydrogenase subunit 5-A, mitochondrial
Expression Region	24-163
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.