



Recombinant Rat ETS-related transcription factor Elf-3 (Elf3)

Product Code	CSB-MP681177RA
Abbreviation	Elf3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4V7E1
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	MAATCEISNV FSNYFNAMYS SEDPTLAPAP LTTFGTEDFV LTLNNQHMSPEGPVGCVGPQ TRSQRDRTEP PAVLHLAEKA SWTGERPQFW SKTQVLWDWIS YQVEKNKYDA SSIDFSRCMDGATLCNCAL EELRLVFGPL GDQLHAQLRD LTSSSSDELS WIIELLEKDG MTFQEGLGDS GPFDQGSPFA QELLDDGRQA SPYYGSSYGP GAPSPGSSDF STSGTDTPQS SHSSDSGGSD VDLDLTDSKV FPRDGFDPYK KGEPKHGKRK RGRPRKLSKEYWDCLEGKKS KHAPRGTHLW EFIRDILHP ELNEGLMKWE NRHEGVFKFL RSEAVAQLWG QKKKNSNMTY EKLSRAMRY YKREILERV DGRRLVYKFGK NSSGWKEEEV GESQN
Source	Mammalian cell
Target Names	Elf3
Protein Names	Recommended name: ETS-related transcription factor Elf-3 Alternative name(s): E74-like factor 3
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.