



Recombinant Rat Coiled-coil domain-containing protein 63 (Ccdc63)

Product Code	CSB-EP681188RA
Abbreviation	Ccdc63
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4V8F7
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	MKKHRRKVSE PLPELSEKAK EQLAEAEELRK LRQQFRKMVD SRKSSNFRNQ RMIAEQYKEI ENLKAQAET TLLSLVKSP KNLDLNQKNF TELRFLQTK GDYEALISSM KVLLGELDDK IVQMERKITN QRQIFLKTQE ANNPRKLQKQ IHILETRLNL VTVHFDTMLT SNAQLRKDIE DLLFEKAAAYD HVYQQQLQRRLL QTQKKTMNVA IEQSTQAYEQ RVEAMARMAA MKDRQQKDIS QYNLEIRELE RLYDHENKLNK SFLLVKLNDR TEFEDQAKKQ EAVKIKKHGK KRKGESFESY EVAHLRLLKL AENGLNQLT EDFLAKEEKN FARFTYVTEL NNDMETMHRK TQRIQDDIIN LRSQQQSSHE GTRNILKQME EKLKRTTEDT DIYENRYREM SKTLEYLKNS VEKLFKINC NATEILGKLG ETGKITDSNL QQYFAIIEKK TNDLLLLESF RRLQEAEGPD VDVPPQFVNP FWGGSALLKP PEPLRVVPPV FGADSFSDKV EEVPHRALSW VQWNTLWTTV ASGSWSWRTT FRRNEEKNHR RHCWRRGTR
Source	E.coli
Target Names	Ccdc63
Protein Names	Recommended name: Coiled-coil domain-containing protein 63
Expression Region	1-559
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.