



Recombinant Human Zinc finger MYND domain-containing protein 17 (ZMYND17)

Product Code	CSB-BP681195HU
Abbreviation	ZMYND17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q4VC12
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAPRSRRRRH KKPPSSVAPI IMAPTTIVTP VPLTPSKPGP SIDTLGFFSL DDNVPGLSQL ILQKLNKMSY EYKLVVDGG TPVSGFGFRC PQEMFQRMED TFRFCAHCRA LPSGLSDSKV LRHCKRCRNV YYCGPECQKS DWPAHRRVCQ ELRLVAVDRL MEWLLVTGDF VLPSGPWPWP PEAVQDWDSW FSMKGLHLDA TLDAVLVSHA VTTLWASVGR PRPDPDLVQG SLKRLLDVLR SRPLTLGLGL RALGIDVRRRT GGSTVHVVGA SHVETFLTRP GDYDELGYMF PGHLGLRVVM VGVDVATGFS QSTSTSPLEP GTIQLSAHRG LYHDFWEEQV ETGQTHHPDL VAAFHPGFHS SPDLMEAWLP TLLLLRDYKI PTLITVYSHQ ELVSSLQILV ELDTHITAFG SNPFMSLKPE QVYSSPNKQP VYCSAYYIMF LGSSCQLDNR QLEEKVDGGI
Source	Baculovirus
Target Names	MSS51
Protein Names	Recommended name: Zinc finger MYND domain-containing protein 17
Expression Region	1-460
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.