



Recombinant Danio rerio Probable thiopurine S-methyltransferase (tpmt)

Product Code	CSB-MP681389DIL
Abbreviation	tpmt
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q504A5
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MSAQADRVMD LSEWENRWQE GRTGFHRSDV HNLLKANVDK LICGRREVRV FFPLCGKAVD MKWLADMGT VVGVEFSEKG IKEFSQEQNL EYNEEAVADI PGAKLFKSTD GKISIQCDL YKFSSAVAGH FGGIWDRGAL VAINPCDRQK YASLLVSLMS SDCRYLLDTL EYNPELYKGP PFFVSEDDIK TVFGGSCNID LLQSVDFGEE KHRSWGLDSL TEKLYLLTTK TQ
Source	Mammalian cell
Target Names	tpmt
Protein Names	Recommended name: Probable thiopurine S-methyltransferase Short name= Thiopurine methyltransferase EC= 2.1.1.67
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.