



Recombinant Human Uncharacterized protein C12orf68 (C12orf68)

Product Code	CSB-BP681476HU
Abbreviation	C12orf68
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q52MB2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEDGLLEIMT KDGGDMPAPL EVSTVPAVGD VISGEYNGGM KELMEHLKAQ LQALFEDVRA MRGALDEQAS HIQVLSDDVC ANQRAIVSMC QIMTTAPRQG GLGVVGGKGS FQSDPQEPET PSPGIGDSSL LGRDPEDEEE EEEEKEMPSP ATPSSHCHERP ESPCAGLLGG DGPLVEPLDM PDITLLQLEG EASL
Source	Baculovirus
Target Names	CCDC184
Protein Names	Recommended name: Uncharacterized protein C12orf68
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.