



Recombinant Human Leucine-rich repeat-containing protein 23 (LRRC23)

Product Code	CSB-MP681504HU
Abbreviation	LRRC23
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53EV4
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSDEDDLEDS EPDQDDSEKE EDEKETEEGE DYRKEGEEFP EEWLPTPLTE DMMKEGLSLL CKTGNGLAHA YVKLEVKERD LTDIYLLRSY IHLRYVDISE NHLTDLSP LN YLTHLLWLKA DGNRLRSAQM NELPYLQIAS FAYNQITDTE GISHPRLETL NLKGNSIHMV TGLDPEKLIS LHTVELRGNQ LESTLGINLP KLKNLYLAQN MLKKVEGLED LSNLTTLHLR DNQIDTSLGF SREMKSLQYL NLRGNMVANL GELAKLRDLP KLRALVLLDN PCTDETSYRQ EALVQMPYLE RLDKEFYEEE ERAEADVIRQ RLKEEKEQEP EPQRDLEPEQ SLI
Source	Mammalian cell
Target Names	LRRC23
Protein Names	Recommended name: Leucine-rich repeat-containing protein 23 Alternative name(s): Leucine-rich protein B7
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.