



# Recombinant Human U4/U6.U5 tri-snRNP-associated protein 2 (USP39)

<b>Product Code</b>	CSB-MP681507HU
<b>Abbreviation</b>	USP39
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q53GS9
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSGRSKRESR GSTRGKRESE SRGSSGRVVKR ERDREREPEA ASSRGSPVRV KREFEPASAR EAPASVVPFV RVKREREVDE DSEPEREVRA KNGRVDSEDR RSRHCPYLDT INRSVLDLDFD EKLCSSISLSH INAYACLVCG KYFQGRGLKS HAYIHSVQFS HHVFLNLHTL KFYCLPDNYE IIDSSLEDIT YVLKPTFTKQ QIANLDKQAK LSRAYDGTTY LPGIVGLNNI KANDYANAVL QALSNVPLR NYFLEEDNYK NIKRPPGDIM FLLVQRFGEL MRKLWNPRNF KAHVSPHEML QAVVLCCKT FQITKQGDGV DFLSWFLNAL HSALGGTKKK KKTIVTDVFQ GSMRIFTKKL PHPDLPAEEK EQLLHNDEYQ ETMVESTFMY LTLDLPTAPL YKDEKEQLII PQVPLFNILA KFNNGITEKEY KTYKENFLKR FQLTKLPPYL IFCIKRFTKN NFFVEKNPTI VNFPIITNVDL REYLSEEVQA VHKNTTYDLI ANIVHDGKPS EGSYRIHVLH HGTGKQWYELQ DLQVTDILPQ MITLSEAYIQ IWKRRDNDDET NQQGA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	USP39
<b>Protein Names</b>	Recommended name: U4/U6.U5 tri-snRNP-associated protein 2 Alternative name(s): Inactive ubiquitin-specific peptidase 39 SAD1 homolog U4/U6.U5 tri-snRNP-associated 65 kDa protein Short name= 65K
<b>Expression Region</b>	1-565
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.