



Recombinant *Oryza sativa* subsp. *japonica* Xylanase inhibitor protein 2 (Chib3H-h)

Product Code	CSB-BP681511OFG
Abbreviation	Chib3H-h
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53NL5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	DDP GLAVYWGRHK EEGSLREACD TGRYTTVIIT FYNAFGHGGRY SLDISGHPLA AVGADIKHCQ SRGITVLLSI GGQGGAYSLP TNASAADVAD NLWNAYLGGH RAGVARPFGD DAAVDGIDFF IDQGGADHYD DLARRLDGYN KYYRGRVGVLTATTRCSYP DHRLEKALAT GVFARIHVRLM FGDEQCTMSP RYSWEKWAAA FPGSKVYIGL VASPEQDSAW MFQKDLYYEM LQFVRSPLNY GGLAIYDRYF DKKANYTGEG
Source	Baculovirus
Target Names	Chib3H-h
Protein Names	Recommended name: Xylanase inhibitor protein 2 Alternative name(s): Class III chitinase homolog h
Expression Region	28-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.