



# Recombinant Human Rho GTPase-activating protein 15 (ARHGAP15)

<b>Product Code</b>	CSB-MP681513HU
<b>Abbreviation</b>	ARHGAP15
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q53QZ3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MQKSTNSDTS VETLNSTRQG TGAVQMRIKN ANSHHDRLSQ SKSMILTDVG KVTEPISRHR RNHSQHILKD VIPPLEQLMV EKEGYLQKAK IADGGKCLRK NWSTSWIVLS SRRIEFYKES KQQALS NMKT GHKPEVDLC GAHIEWAKEK SSRKNVFQIT TVSGNEFLLQ SDIDFIILDW FHAIKNAIDR LPKDSSCPSR NLELFKIQRS SSTELLSHYD SDIKEQKPEH RKSLMFRLHH SASDTS DKNR VKSRLKKFIT RRPSLKTLE KGLIKDQIFG SHLHKVCERE NSTVPWFVKQ CIEAVEKRGL DVDGIYRVSG NLATIQKLRF IVNQEEKLNL DDSQWEDIHV VTGALKMFFR ELPEPLFPYS FFEQFVEAIK KQDNNTIEA VKSLVQKLP PNRDTMKVLF GHLTKIVAKA SKNLMSTQSL GIVFGPTLLR AENETGNMAI HMYVQNQIAE LMLSEYSKIF GSEED
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ARHGAP15
<b>Protein Names</b>	Recommended name: Rho GTPase-activating protein 15 Alternative name(s): ArhGAP15 Rho-type GTPase-activating protein 15
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.