



Recombinant Human Uncharacterized protein C2orf53 (C2orf53)

Product Code	CSB-EP681514HU-B
Abbreviation	C2orf53
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53SZ7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MLPQNKDQVL PQTSVLPGRP TWGFSQLVDS SPHNLQPLSP HQGLPPSQPP FSSTQSRPSP SPPASPSPG FQFGSCDSNS DFAPHPYSPS LPSSPTFFHQ NYLSLPRPRA SSPSNHWLYP SPPLTPSFSP SQPQNSSLPH SPCQSPSHPE ELHSSTLTSP GPSPPSHRLH SNRQTRWRHQ YRDTGSGSPG VVERCVPSEK DPAQFRDPGA LAQALVVQLG HRRIAHDRL LLLQHLWLGR TGQAPVVEYP ICLVCLRP RS PSCPLPRYRT GPRLLAFPQL LPCVQGQESG PLRIGIGFGL RLPQGQARAL HLLPEKRPKE AGPQGKATQA CGHQLPASQP PAAQARADPV PGTPSQTRSF RSAGLQSPNS PRCFSGPPPR APKQVTTSLK PRPCPGPKRP VSLELILQKS SV
Source	E.coli
Target Names	PRR30
Protein Names	Recommended name: Uncharacterized protein C2orf53
Expression Region	1-412
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.