



Recombinant Dictyostelium discoideum Mitotic spindle assembly checkpoint protein MAD2B (mad2l2)

Product Code	CSB-MP681717DKK
Abbreviation	mad2l2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54S00
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MLDEHFCDCI GEFLETSFHC ILYIRGVYPS CLFSKSIKYD IPVPISRSDL LTRYISNSID SLKPHFLKDT IEKISLTILN KYDKPIEKFI FEISSLNNTT NNKNKNNNNN NNNSDDDET F NYYINIQQQN YTSNNNNNNN NNINNNNNKP NLLQLEASFK AYIIKILMTT DSFYSEQNFY KTQQSSSSSS SLKTSGGGDG NGFISNSDGV IHFNEIDNET YQIKSNQKN KKEDNDDNNN GDKDLKFTIH VHTKPTTSNI GLNTNLSNVI LNINPLESPR RTNTTINQIT PVWISTTEKD SEIENSYIIP LNSNITDGKT TIRTFTEQSI TK
Source	Mammalian cell
Target Names	mad2l2
Protein Names	Recommended name: Mitotic spindle assembly checkpoint protein MAD2B Alternative name(s): Mitotic arrest deficient 2-like protein 2 Short name= MAD2-like protein 2
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.