



Recombinant Dictyostelium discoideum 40S ribosomal protein S6 (rps6)

Product Code	CSB-BP681751DKK
Abbreviation	rps6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54UU3
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MKLNIAYPVN GSQKKINIED KNKVRCFMEK RIGQEVEATT LGDEFKGYVF KITGGNDTEG FPMMQGVGAP TRVRLLLDGR SGCFKPSRDG ERKRKSVRGC IVAEDIASLQ LIIVKKGDAE IPGLTDVSFP ASKGPKRASN IRKLFKLSKD EDVRSFVIRR ELPATDKRKA KSKAPKIQRL VTPTTVARRK ALRAKKAARH TKASNEAAEY AKLVAQRQQA KAKRSVKVSS KKVVAK
Source	Baculovirus
Target Names	rps6
Protein Names	Recommended name: 40S ribosomal protein S6
Expression Region	1-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.