



Recombinant Synechocystis sp. Light-dependent protochlorophyllide reductase (por)

Product Code	CSB-EP682427SSQ-B
Abbreviation	por
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q59987
Product Type	Recombinant Protein
Immunogen Species	Synechocystis sp. (strain PCC 6803 / Kazusa)
Purity	>85% (SDS-PAGE)
Sequence	MEQPMKPTVI ITGASSGVGL YGAKALIDKG WHVIMACRNL DKTQKVADEL GFPKDSYTII KLDLGYLDSV RRFVAQFREL GRPLKALVCN AAVYFPLLDE PLWSADDYEL SVATNHLGHF LLCNLLLEDL KACPDADKRL IILGTVTANS KELGGKIPIP APPDLGNFEG FEAGFKKPIA MINNKKFKSG KAYKDSKLCN MLTTRELHRR FHQETGIVFN SLYPGCVADT PLFRNHYSLF RTIFPWFQKN VTKGYVSQEL AGERVAMVVA DDKFKDSGVH WSWG NRQQAG REAFVQELSE QGSDAQKAQR MWDLSEKLVG LV
Source	E.coli
Target Names	por
Protein Names	Recommended name: Light-dependent protochlorophyllide reductase Short name= PCR EC= 1.3.1.33 Alternative name(s): NADPH-protochlorophyllide oxidoreductase Short name= LPOR Short name= POR
Expression Region	1-322
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.