



Recombinant Candida albicans Mitochondrial distribution and morphology protein 12 (MDM12)

Product Code	CSB-EP682442CZD
Abbreviation	MDM12
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q59S52
Product Type	Recombinant Protein
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSFDINWNQL TIDDTINQSI KEFLDQQFKN ISLPSFISNL AVTDFNLGEI PPEVTIRHIG DPFEFYEDE>NNLGMNEAN NDSKDEHLKN HGDGINKDSG YNSQNLDDDED EDEDEDDDED EDEEEDEDD YDDHDLGTIN EGISLLNFNE NSTTPSANSF AGSAAPPLPP PLNPSRDSFH SILHPYGVNS IIGATGAGSE TPTNILNQNY LSSRVLPKIS VKQKQPHHDD NDIQLIVEIN YKGDMMHINLL VNLLVNYPSP NFISLPIKLH ITDIVIHSIA TIAYLKKS VF LSFLCDVDDT FPDFDSNVQT PTSTTGGNFV DYYSNDATIN KERIDIVKKI KIESEIGEVE NNILRNVGKV EKFLVEQLRN ILRDEIAWPS WICIDMND DD DEEEEEEESE DNDGGNSDLN DNDGKHGDGR TDETEAGE
Source	E.coli
Target Names	MDM12
Protein Names	Recommended name: Mitochondrial distribution and morphology protein 12 Alternative name(s): Mitochondrial inheritance component MDM12
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.