



Recombinant Emericella nidulans Pre-mRNA-splicing factor slu7 (slu7)

Product Code	CSB-BP682522EKF
Abbreviation	slu7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5B3U2
Product Type	Recombinant Protein
Immunogen Species	Emericella nidulans (strain FGSC A4 / ATCC 38163 / CBS 112.46 / NRRL 194 / M139) (Aspergillus nidulans)
Purity	≥85% (SDS-PAGE)
Sequence	MSRKPADVAS KERNEYIPAF ISKKPFYIDD DDTANDYLEH QRLHKQTTDQ SKWYER GKRA GPAATKYRKG ACENCGAMTH KAKECLSRPR KHGAKWTGKD IQADEVIQNV DLGWDAKRDR WNGYDAAEYR QVVEEYEELE RLKRQAKLTK GETQTTNDGL DDEAPEQEAR YAEESDMGRQ QSTATRNLRI REDTAKYLLN LDLDSAKYDP KTRRMVDMGA AEDQAAALVA EENFVRSSGD AAEFERAQR Y AWEAQERG TQ KIH LQANPTS GEITRKKELA ESEAKRDAHR KALLEKYGGE QHLKHTPLLE TMVVENERFV EYDETGAIKG APKKATKSKY PEDILTNNHK SVWGSWWHNF QWGYACCFST VKNSYCTGEE GKRAFEEARN MLLLPGDETE QPSLAVESAS RQEEPSAESH NQQRDSKKRT LMEVQSGITE EELESYKRSR LAADDPMAAF IEKDDS
Source	Baculovirus
Target Names	slu7
Protein Names	Recommended name: Pre-mRNA-splicing factor slu7 Alternative name(s): Splicing factor sluA
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.