



Recombinant Emericella nidulans Flap endonuclease 1 (fen1)

Product Code	CSB-EP682535EKF
Abbreviation	fen1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5B9L6
Product Type	Recombinant Protein
Immunogen Species	Emericella nidulans (strain FGSC A4 / ATCC 38163 / CBS 112.46 / NRRL 194 / M139) (Aspergillus nidulans)
Purity	>85% (SDS-PAGE)
Sequence	MGIKHLYQVI AENAPDAIKA GDIKNHFGRK VAIDASMSIY SFLIAVRSEG QQLMSDTGET TSHLMGMFYR TLRMVDNGIK PLYVFDGAPP KLKSGELAKR SARKHEATEA HEEAKETGTA EDVEKFSRRT VRVTREHNAE CKKLLKLMGI PYIDAPTEAE AQCAVLARAG KVYAAASEDM DTLCFEAPIL LRHLTFSEQR KEPIQEIHLN RALEGLGMDR NQFIDLCILL GCDYLEPIPK VGPNTALKLI RDHGTLEKVV ESIEKDPKQK YVIPESWPYQ DARELFLNPD VRDANDPECD FKWEAPDIEG LVDFLVKDKG FNEDRVRNGA ARLQKNLKTA QQRLEGFFK PVAKTDAEKA TLKRKHDEKI QEQQKRRKKEE AKAKKEAKAR PRAAV
Source	E.coli
Target Names	fen1
Protein Names	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.