



Recombinant Danio rerio PHD finger protein 23A (phf23a)

Product Code	CSB-BP682558DIL
Abbreviation	phf23a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BJ10
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MLGIMDHHQD TVRKCKSEAL PPERRKRTVE DFNKFCFVL TYAGYIPPQK EESSWSPSSS PCTHDLSELS GEGSVKDSWT DSHSDLNNIH NLVYKAETDS SSSREFSHLP SDNSLDKMTL KDSLNVHHSK AERKKVKKLD RLSLGGPRKS ISEARAHKQS KAALKKIKTS IKAERHFTSS SPLNEGFEKE ELTEQAIHMH EAGLKLESNQ ETDLSSCETD TLVTDEDIMV ESGDSDWDLI TCYCGKPFAG RPMIECEECS IWWHLSCAKI KKSNDPDIY CYRCLDSRGS TVKRDH
Source	Baculovirus
Target Names	phf23a
Protein Names	Recommended name: PHD finger protein 23A
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.