



Recombinant Danio rerio Methyltransferase-like protein 21A (mettl21a)

Product Code	CSB-YP682579DIL
Abbreviation	mettl21a
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BLD8
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MALVPYDENV LPALSKLHQS SAEFTLANHR IRLSQDWKRL GVAAVVWDAA VVLCMFLEMV KVDLKGKRVV ELGAGTGLVG IVAALLGANV TITDREPALE FLTANVHENV PQGRQKAVQV SELTWGENLD LYPQGGYDLI LGADIVYLEE TFPALLQTLE HLSSGDTVVL LSCRIRYERD ERFLTELRRR FSVQEVHYDS QRDIHVYRAV KNKSNTLE
Source	Yeast
Target Names	mettl21a
Protein Names	Recommended name: Methyltransferase-like protein 21A EC= 2.1.1.-
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.