



Recombinant Bovine Cell cycle checkpoint control protein RAD9B (RAD9B)

Product Code	CSB-EP682705BO-B
Abbreviation	RAD9B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5E9X8
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	≥85% (SDS-PAGE)
Sequence	MKAHLLLFMR IYIHRCP LN EVAPPGQVVE VFGKAIQALS RVSDDELWLDP SEKGLALRSV NSCRSAYGCV LFSPVFFQHY QWSASVKMND TDIILNLNCR LGMKSILPIF RCLNSLEKNV EKCKIFTRSD KCKVVIQFFC RHGIKKIHNV CFQGSRPLQV IFQKNMCANT LVIQPRVLA E AIVLFTSSQE EVTLAVTPLK VCIKSSNEES MDLTDSVYSE MFVGPDEFDF FQIGIDTEIT FCFKELKGV L TFSEAIHAPI AIHFDFPGKP MALSIDDMLL EANFILATLA DEPSRASSLQ TLYLSQKQRR SEPIHSNSKA GKNITSKVPE YISRKVEPKR LYSNETPTNI STLENCGSPL MKRANKDITE VPESDGNLSE VPESVSDTE DVPGSPCLKK FSCMFFGAAS SDQQEPFSLP FQSLATASGS EEDMNNNGSFS TF
Source	E.coli
Target Names	RAD9B
Protein Names	Recommended name: Cell cycle checkpoint control protein RAD9B Alternative name(s): DNA repair exonuclease rad9 homolog B
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.