



# Recombinant Human Arylsulfatase J (ARSJ)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-MP682915HU                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Abbreviation</b>      | ARSJ                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Uniprot No.</b>       | Q5FYB0                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Immunogen Species</b> | Homo sapiens (Human)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Purity</b>            | ≥85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Sequence</b>          | L EEEEEGALLA QAGEKLEPST TSTSQPHLIF ILADDQGFRD VGYHGSEIKT<br>PTLDKLA AEG VKLENYYVQP ICTPSRSQFI TGKYQIHTGL QHSIIRPTQP<br>NCLPLDNATL PQKLKEVGYS THMVGKWHLG FYRKECMPTR RGFDTFFGSL<br>LGSGDY YTHY KCDSPGMCGY DLYENDNAAW DYDNGIYSTQ MYTQRVQQIL<br>ASHNPTKPIF LYIAYQAVHS PLQAPGRYFE HYRSIININR RRYAAMLSC L<br>DEAINNVTLA LKTYGFYNNS IIIYSSDNGG QPTAGGSNWP LRSGSKGTYWE<br>GGIRAVGFVH SPLLNKKGTV CKELVHITDW YPTLISLAEG QIDEDIQLDG<br>YDIWETISEG LRSPRVDILH NIDPIYTKAK NGSWAAGYGI WNTAIQSAIR<br>VQHWKLLTGN PGYSDWVPPQ SFSNLGPNRW HNERITLSTG KSVWLFNITA<br>DPYERVDLSN RYPGIVKKLL RRLSQFNKTA VPVRYPPKDP RSNPRLNGGV<br>WGPWYKEETK KKKPSKNQAE KKQKKS KKKK KKQQKAVSGS TCHSGVT CG |
| <b>Source</b>            | Mammalian cell                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| <b>Target Names</b>      | ARSJ                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| <b>Protein Names</b>     | Recommended name: Arylsulfatase J Short name= ASJ EC= 3.1.6.-                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Expression Region</b> | 50-599                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Protein Length</b>    | Full Length of Mature Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                                                                                                       |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                                                                                                |