




Recombinant Staphylococcus aureus

Staphylococcal secretory antigen ssaA1 (ssaA1)

Product Code	CSB-YP683042FLB
Relevance	Not known; immunogenic protein.
Abbreviation	Recombinant Staphylococcus aureus ssaA1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5HCY4
Product Type	Recombinant Protein
Immunogen Species	Staphylococcus aureus (strain COL)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	AEQNNNGYNSNDAQSYSYTYTIDAQGNYHYTWTGNWNPSQLTQNNTYYYNN YNTYSYNNASYNNNYNHSYQYNNYTNNSTATNNYYTGGSYGASTTSSNNV HVTTTAAPSSNGRSISNGYASGSNLVTSGQCTYYVFDRVGGKIGSTWGNASN WANAAASSGYTVNNTPKVGAIMQTTQGYGHHVAYVEGVNSNGSVRVSEMNY GHGAGVVTSTRTISANQAGSYNFIH
Research Area	Others
Source	Yeast
Target Names	ssaA1
Protein Names	Recommended name: Staphylococcal secretory antigen ssaA1
Expression Region	27-255aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	27.0kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.