



Recombinant Pyrococcus kodakaraensis Cell division control protein 6 homolog (cdc6)

Product Code	CSB-BP683324FHY
Abbreviation	cdc6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JET2
Product Type	Recombinant Protein
Immunogen Species	Thermococcus kodakarensis (strain ATCC BAA-918 / JCM 12380 / KOD1) (Pyrococcus kodakaraensis (strain KOD1))
Purity	≥85% (SDS-PAGE)
Sequence	MDDDYLSSIF EKYLHARKIF KNKEVLRHSY TPKELPHRHK QIDDLAHLV PVL RGETPSN VFVYGKTGTG KTVTVKFVTE ELKKISQKYN IPVEVIYINC EIIDTHYRVL ARIVNHFKEE SGIEVPLVGW PTDEVYAKLK EVIDAKERFV IIVLDEIDKL IKKSGDDILY SLTRINSELS KAKVSIIGIS NDLKFKEYLD ARVLSSLSEE EVVFPYDAN QLRDILMQRA KEAFYEGVLD DAVVPLCAAL AAREHGDARR ALDLLRVAGE IAEREGASKV TERHVWKAQE KIEQDTMEEV IKTLPLHSKV LLYAIVMLDE NGELPANTGE VYSIYKMLCD SLDVEPLTQR RVSDLINELD MLGIINAKVV SKGRYGRTE IRLNVTPYMV KNIYRHDEQV RPLLTLSLSR QRRLF
Source	Baculovirus
Target Names	cdc6
Protein Names	Recommended name: Cell division control protein 6 homolog Short name= CDC6 homolog
Expression Region	1-415
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.