



Recombinant *Oryza sativa* subsp. *japonica* Auxin-responsive protein IAA3 (IAA3)

Product Code	CSB-EP683701OFG
Abbreviation	IAA3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5NB25
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSPPLELDYI GLSPPPPPPS SSSAAAARAD DVDLKGTELRL LGLPGSESPD RRPAAIAAAA ATATTLELLP AKGAKRVFPD EAALTPPTAA AGKGKAAREG EEVGAEEDK KVAAPPQPAA KAQVVGWPPV RSYRKNTMAT NQIKSNKEDV DAKQGQGFY VKVSMDGAPY LRKVDLKYK NYKDMSLGLK KMFIFSTGK EGAENQKDGE YVLTIEDKDG DWMLVGDVPW EMFTDSCRRL RIMKGSDAIG LAPRAGEKSK NRN
Source	E.coli
Target Names	IAA3
Protein Names	Recommended name: Auxin-responsive protein IAA3 Alternative name(s): Indoleacetic acid-induced protein 3
Expression Region	1-263
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.