



# Recombinant Danio rerio Nuclear receptor subfamily 6 group A member 1-B (nr6a1b)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-EP684149DIL  |
| <b>Abbreviation</b>      | nr6a1b   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q4V8R7   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Danio rerio (Zebrafish) (Brachydanio rerio)  |
| <b>Purity</b>            | >85% (SDS-PAGE)  |
| <b>Sequence</b>          | MEVENRPNTY DFSKRDKLLS NRFYEENSLQ DETDDGGERW CLICGDRASG<br>LHYGIISCEG CKGFFKRSIC NKRIYRCNRD KNCQMSRKQR NRCQYCRLQK<br>CLQMGMNRKA IREDGMPGGR NKMIGPVHIS LEEIERLMSG QEFKEGSDLS<br>DSWSHGYSNH SSPGNSLSEG GQSLSFSSSR SVSCRDECIS PQLTHTFLMC<br>KYPLPPPTGS ILKTQTHLT GQILADEDLT PLTPMLIED GYSVTQSELL<br>ALLCGIADEL LFRQIVWLKR LPFFTDLSIK DCTRLLGSSW HQLILLSSIT<br>VHSAQILGEL ANVTHHYTPS SHTLQRFGE AMEVMESLNF LFRKFHQLNI<br>SNEEYSCLKT ITLLNQETTG LCNTSMLKQL SERYWTLCRE LTERLHPQRP<br>KRFSDIITCL TEIRHTSGKM MSIPLEQLPL LFKAVLYSCT TNQNPWLPKS<br>STSRT |
| <b>Source</b>            | E.coli   |
| <b>Target Names</b>      | nr6a1b   |
| <b>Protein Names</b>     | Recommended name: Nuclear receptor subfamily 6 group A member 1-B<br>Alternative name(s): Germ cell nuclear factor B Short name= GCNF-B  |
| <b>Expression Region</b> | 1-455  |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | full length protein  |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life  |



of lyophilized form is 12 months at -20°C/-80°C.