



# Recombinant Plasmodium berghei Flap endonuclease 1 (FEN1)

<b>Product Code</b>	CSB-BP684246PLH
<b>Abbreviation</b>	FEN1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q4Z015
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Plasmodium berghei (strain Anka)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	<p>           MGIKGLTKFI ADTAPNAIKE IKIENLMGRV VAIDASMSLY QFIIAIRDGD            QYGNLMNESG ETTSHISGLM SRTIKLMENG LKPIYVFDGA PPELKGSELE            KRGEKRQKAE ELLIKAKEEN LEEIKKQSGR TVRVTKKQNE EAKKLLTLMG            IPVIESPCEA EAQCAFLTKY EMAHATATED ADALVFGTKI LIRNLNANAS            SNKNKNKNSS KRGYILTEIN LEQVLKGLKL TMDEFIDFCI LCGCDYCDTI            KGIGSKTAYN LIKEYNCIEN IIKNIDQNKY QVPANFKYVE ARQSFIPKV            LEKSEVKIDW CEPKIEELKK FLIKEHNFNE VRVTNYITRL LKARKVTTQR            RLDTFFTTCT KKSTKLIIEE NQKEVLKTKS KGKKREINND NSTKLNAKKK            KTNVKDEKKN NEKVDEPKNK PNENLVKDED DQDDYDQNIIF DERTNSESGN            IKNENVKEDI SSSDITMDIP KCINDIVC         </p>
<b>Source</b>	Baculovirus
<b>Target Names</b>	FEN1
<b>Protein Names</b>	Recommended name: Flap endonuclease 1 Short name= FEN-1 EC= 3.1.-.- Alternative name(s): Flap structure-specific endonuclease 1
<b>Expression Region</b>	1-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.