



Recombinant RepFIB replication protein A (repA)

Product Code	CSB-EP684409ENL-B
Abbreviation	repA
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q52219
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli
Purity	>85% (SDS-PAGE)
Sequence	MDKSSGELVT LTPNNNTVQ PVALMRLGVF VPTLKSLKNS KKNTLSRTDA TEELTRLSLA RAEGFDKVEI TGPRLDMDND FKTWVGIIHS FARHNVIGDK VELPFVEFAK LCGIPSSQSS RRLRERISPS LKRIAGTVIS FSRTDEKHTR EYITHLVQSA YYDTERDIVQ LQADPRLFEL YQFDRKVLLQ LKAINALKRR ESQAQALYTFI ESLPRDPAPI SLARLRARLN LKSPVFSQSQ TVRRAMEQLR EIGYLDYTEI QRGRTKLFCL HYRRPRLKAP NDESKENPLP PSPAEKVSPE MAEKLALLEK LGITLDDLEK LFKSR
Source	E.coli
Target Names	repA
Protein Names	Recommended name: RepFIB replication protein A
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>