



Recombinant Mouse Cyclin-L1 (Ccnl1)

Product Code	CSB-MP684434MO
Abbreviation	Ccnl1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q52KE7
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MASGPHPTST AAAAAAAAAA ASSAAPSAGG SSSGTTTTTT TTTGGILIGD RLYSEVSLTI DHSLIPEERL SPTPSMQDGL DLPSETDLRI LGCELIQAAG ILLRLPQVAM ATGQVLFHRF FYSKSFVKHS FEIVAMACIN LASKIEEAPR RIRDVINVFH HLRQLRGKRT PSPLILDQNY INTKNQVIKA ERRVLKELGF CVHVKHPHKI IVMYLQVLEC ERNQLTVQTA WNYMNDLRT NVFVRFQPET IACACIYLA RALQIPLPTR PHWFLFGTT EEEIQEICIE TLRLYTRKKP NYELLEKEVE KRKVALQEAK LKAKGLNLDG TPALSTLGGF SPASKPSSPR EVKAEKSPV SINVKTVKKE PEDRQQASKS PYNGVRKDSK RSRTSRSASR SRSRTRSRSR SHSPRRHYNN RRSRSGTYSS RSRSRSRSHS ESPRRHHNHG SPHLKAKHTR EDLKSSNRHG HKRKKRSRSR QSKTRDHSV TTKHRHERGH HRDRRERSRS FERSHKGKHH GGSRSRGHRH RR
Source	Mammalian cell
Target Names	Ccnl1
Protein Names	Recommended name: Cyclin-L1 Short name= Cyclin-L Alternative name(s): Cyclin Ania-6a
Expression Region	1-532
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.