



Recombinant Human Protein FAM98B (FAM98B)

Product Code	CSB-MP684437HU
Abbreviation	FAM98B
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q52LJ0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRGPEPGPQP TMEGDVLDLTL EALGYKGPLL EEQALTKAAE GGLSSPEFSE LCIWLGSQIK SLCNLEESIT SAGRDDLESF QLEISGFLKE MACPYSVLIS GDIKDRLLKKK EDCLKLLLFL STELQASQIL QNKKHKNSQL DKNSEVYQEV QAMFDTLGIP KSTTSDIPHM LNQVESKVVD ILSKVQKNHV GKPLLKMDLN SEQAEQLERI NDALSCEYEC RRRMLMKRLD VTVQSFGWSD RAKVKTDDIA RIYQPKRYAL SPKTTITMAH LLAAREDLSK IIRTSSGTSR EKTACAINKV GVSFSTVENE LMISYLMFLQ ILVYFSFMSW
Source	Mammalian cell
Target Names	FAM98B
Protein Names	Recommended name: Protein FAM98B
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>