



# Recombinant Human Kelch-like protein 22 (KLHL22)

<b>Product Code</b>	CSB-MP684469HU
<b>Abbreviation</b>	KLHL22
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q53GT1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	AEEQEFTQL CKLPAQPSHP HCVNNTYRSA QHSQALLRGL LALRDSGILF DVVLVVEGRH IEAHRILLAA SCDYFRGMFA GGLKEMEQEE VLIHGVSINA MCQILHFIYT SELELSLSNV QETLVAACQL QIPEIIHFCC DFLMSWVDEE NILDVYRLAE LFDLSRLTEQ LDTYILKNFV AFSRTDKYRQ LPLEKVYSLL SSNRLEVSCE TEVYEGALLY HYSLEQVQAD QISLHEPPKL LETVRFPLME AEVLQRLHDK LDPSPLRDTV ASALMYHRNE SLQPSLQSPQ TELRSDFQCV VGFGGIHSTP STVLSDQAKY LNPLLGEWKH FTASLAPRMS NQGIAVLNNF VYLIGGDNNV QGFRAESRCW RYDPRHNRWF QIQSLQQEHA DLSVCVVGRY IYAVAGRDYH NDLNAVERYD PATNSWAYVA PLKREVYAHA GATLEGKMYI TCGRRGEDYL KETHCYDPGS NTWHTLADGP VRRAWHGMAT LLNKLYVIGG SNNDAGYRRD VHQVACYSCT SGQWSSVCPL PAGHGEPGIA VLDNRIYVLG GRSHNRGSRT GYVHIYDVEK DCWEEGPQLD NSISGLAACV LTLPRSLLE PPRGTPDRSQ ADPDFASEVM SVSDWEEFDN SSED
<b>Source</b>	Mammalian cell
<b>Target Names</b>	KLHL22
<b>Protein Names</b>	Recommended name: Kelch-like protein 22
<b>Expression Region</b>	2-634
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.