



Recombinant *Oryza sativa* subsp. *japonica* Putative B3 domain-containing protein Os11g0242900 (Os11g0242900, LOC_Os11g13900)

Product Code	CSB-MP684474OFG
Abbreviation	Os11g0242900, LOC_Os11g13900
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53N90
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTVELEKIAG SFFISKGWKT FVHRTGLLSG QYIRFQVLTP SKINVLLFDK KKDSKLPMP SSKKQIKTAP KRSTGITIND MPTSKHASML ISHTSNKETS SDSRTESMTD IPSSSDNSGE TTRSFDDLFC CARNTAVTPD IKNYISIIGQ FLQRSSKFYI VTMNNTFMKQ DRLVEAAGSD SVTMLLHKSS DDRCNLKRGW ATFAATNAIH LHSVCIFHFY KPPNVKITID VL
Source	Mammalian cell
Target Names	Os11g0242900
Protein Names	Recommended name: Putative B3 domain-containing protein Os11g0242900
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.