



Recombinant Human Retroviral-like aspartic protease 1 (ASPRV1)

Product Code	CSB-EP684475HU-B
Abbreviation	ASPRV1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q53RT3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	SMGKGYLKGKIGKVPVRFLVDSGAQVSVVHPNLWEEVTDGDLDTLQPFENV VKVANGAEMKILGVWDTAVSLGKCLKLKAQFLVANASAEAAIIGTDVLQDHNAIL DFEHRCTCLKGKKFRLLPVGGSLEDEFDLE
Source	E.coli
Target Names	ASPRV1
Protein Names	Recommended name: Retroviral-like aspartic protease 1 EC= 3.4.23.- Alternative name(s): Skin-specific retroviral-like aspartic protease Short name= SASPase Short name= Skin aspartic protease TPA-inducible aspartic proteinase-like
Expression Region	191-326aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.