



Recombinant Dictyostelium discoideum Inorganic pyrophosphatase (ppa1)

Product Code	CSB-EP684663DKK-B
Abbreviation	ppa1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q54PV8
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	MKNFFSFSKL TPHLKNIYHR NMTYTTKQVG ETGSLEYRLF FLKDNKPVSS FHDVPLWVVK EKQIVNMLVE IPRGTNAKLE IATKEYMNPI KQDVKDGKLR FVHDKYPFNY GALPQTWESP EHTHPSTGAK GDNDPLDACE IGSGQGVTGE FKQVKVLGVF AMIDAGETDW KILCIDVNDP IASQINSQED IEKHLPKIN EVYTFLRDYK IPDGKGPNQF AFDGKLQSID FSMKIIETE AEWKDLVGGK TKSSLSVVNT TLNDSNTVTA DVAAQKLNK
Source	E.coli
Target Names	ppa1
Protein Names	Recommended name: Inorganic pyrophosphatase EC= 3.6.1.1 Alternative name(s): Pyrophosphate phospho-hydrolase Short name= PPase
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.