



Recombinant Dictyostelium discoideum ATP synthase subunit delta, mitochondrial (atp5D)

Product Code	CSB-YP684886DKK
Abbreviation	atp5D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q55F42
Product Type	Recombinant Protein
Immunogen Species	Dictyostelium discoideum (Slime mold)
Purity	>85% (SDS-PAGE)
Sequence	SDLLTFS LLSPHQTIYK DKKAQLVTLP GAKGIFGVAK NHVPRIAE LK PGVIQINHEN GDLEKFFISG GFAPVNP DAS CYINTIEAVP IDQLDAEEVK NGLARYTQLY NDAQEENAKA VALIGLETYQ QMAFACGVSA
Source	Yeast
Target Names	atp5D
Protein Names	Recommended name: ATP synthase subunit delta, mitochondrial Alternative name(s): F-ATPase delta subunit
Expression Region	34-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.