



Recombinant Danio rerio WD repeat-containing protein 73 (wdr73)

Product Code	CSB-MP684944DIL
Abbreviation	wdr73
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q567G2
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MDISSDESDD WMESLNAYR DLHLFQLEHP TKVIEWMGDK SICVAGYTGA RSEILELLLP LKLYAGDNQG LCAERDFKVK HGGFSEEPVE CLIHIPGTRC LVTSGNSNTL QIWDIGGDDS DVIKMTRVIN LKSPSAKGSK ITPGFTDSPT VLHGSCISDI QLTEIATGRV LHTVGKSSD SVSGLQFVNA FVFLICSTNG ALLLGDTRDP STCHYPLEES NSGLRWSFGL RSDSSQSEPS CCKVARLSSS GHMLLSDLRD LKSPFCQTQI KVPQSASKNH FMNISWAPVL DKCLSVSGFD GTVHIYNTNN WSTESSSPQP IFSHRGHEMS QEAKSSTSPA VVTAHSWHPS RPKTLLSAAT DGSLHVWDWV DKITDR
Source	Mammalian cell
Target Names	wdr73
Protein Names	Recommended name: WD repeat-containing protein 73
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.