



# Recombinant Human Coiled-coil domain-containing protein 93 (CCDC93)

<b>Product Code</b>	CSB-BP684945HU
<b>Abbreviation</b>	CCDC93
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q567U6
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MGLPRGPEGQ GLPEVETRED EEQNVKLTEI LELLVAAGYF RARIKGLSPF DKVVGGMTWC ITTCNFDVDV DLLFQENSTI GQKIALSEKI VSVLPRMKCP HQLPHQIQG MDFIHFPVV QWLVKRAIET KEEMGDYIRS YSVSQFQKTY SLPEDDDFIK RKEKAIKTVV DLSEVYKPRR KYKRHQGAEE LLDEESRIHA TLLEYGRRYG FSRQSKMEKA EDKKTALPAG LSATEKADAH EEDELRAAEE QRIQSLMTKM TAMANEE SRL TASSVGQIVG LCSAEIKQIV SEYAEKQSEL SAEESPEKLG TSQLHRRKVI SLNKQIAQKT KHLEELRASH TSLQARYNEA KKT LTE LKTY SEKLDKEQAA LEKIESKADP SILQNLRALV AMNENLKSQE QEFKAHCREE MTRLQQE IEN LKAERAPRGD EKTLS SGEPP GTLTSAMTHD EDLDRRYNME KEKLYKIRLL QARRNREIAI LHRKIDEVPS RAE LIQYQKR FIELYRQISA VHKETKQFFT LYNTLDDKKV YLEKEISLLN SIHENFSQAM ASPAARDQFL RQMEQIVEGI KQSRMKMEKK KQENKMRRDQ LNDQYLELLE KQRLYFKTVK EFKEEGRKNE MLLSKVKAKA S
<b>Source</b>	Baculovirus
<b>Target Names</b>	CCDC93
<b>Protein Names</b>	Recommended name: Coiled-coil domain-containing protein 93
<b>Expression Region</b>	1-631
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.