



Recombinant Human Coiled-coil domain-containing protein 93 (CCDC93)

Product Code	CSB-MP684945HU
Abbreviation	CCDC93
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q567U6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MGLPRGPEGQ GLPEVETRED EEQNVKLTEI LELLVAAGYF RARIKGLSPF DKVVGGMTWC ITTCNFDVDV DLLFQENSTI GQKIALSEKI VSVLPRMKCP HQLPHQIQG MDFIHFPVV QWLVKRAIET KEEMGDYIRS YSVSQFQKTY SLPEDDDFIK RKEKAIKTVV DLSEVYKPRR KYKRHQGAEE LLDEESRIHA TLLEYGRRYG FSRQSKMEKA EDKKTALPAG LSATEKADAH EEDELRAAEE QRIQSLMTKM TAMANEE SRL TASSVGQIVG LCSAEIKQIV SEYAEKQSEL SAEESPEKLG TSQLHRRKVI SLNKQIAQKT KHLEELRASH TSLQARYNEA KKTLTELKTY SEKLDKEQAA LEKIESKADP SILQNLRALV AMNENLKSQE QEFKAHCREE MTRLQQE IEN LKAERAPRGD EKTLSGGEPP GTLTSAMTHD EDLDRRYNME KEKLYKIRLL QARRNREIAI LHRKIDEVPS RAE LIQYQKR FIELYRQISA VHKETKQFFT LYNTLDDKKV YLEKEISLLN SIHENFSQAM ASPAARDQFL RQMEQIVEGI KQSRMKMEKK KQENKMRRDQ LNDQYLELLE KQRLYFKTVK EFKEEGRKNE MLLSKVKAKA S
Source	Mammalian cell
Target Names	CCDC93
Protein Names	Recommended name: Coiled-coil domain-containing protein 93
Expression Region	1-631
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.