



Recombinant Arabidopsis thaliana FBD-associated F-box protein At5g18780 (At5g18780)

Product Code	CSB-BP684972DOA
Abbreviation	At5g18780
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q56YH2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	<pre>MQGRVRGSGE DRISILPEPL LCHILSFLRT KDSVRTSVLS SRWRDLWLWV PRLDLDKSD F SDDNPSASFV DKFLNFRGES YLRGFKLNTD HDVYDISTLD ACLMRLDKCK IQHFEIENCF GFCILLMPLI IPMCHTLVSL KLSFVILSKF GSLSLPCLEI MHFEKVI FPS DKSAEVL IAC SPVLRDLRIS QSGDDAVEVL RVCSASLKS F TLKRTDHDYV GNGEYTVVID TPRLEYLNLK DYQCQGFKIV SMSEYVRVHV DVVFEVIGGT VLSKRNI CD FLSCVSNVRY MTISRRSLEF IYRHLELKPR FKFHDLARLR ATMFSNSSPE MLPVILETCP NLKHLTLELV HDSLVT EGTS GLLTVLPRCL ISSLASVDIE SPITDKATEL KLVSYLLENS TTLKKLV LRL NQSCRDKYEP GLLKQVLQSP RCSSLCQLVI L</pre>
Source	Baculovirus
Target Names	At5g18780
Protein Names	Recommended name: FBD-associated F-box protein At5g18780
Expression Region	1-441
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.