



Recombinant Danio rerio Homeobox protein Hox-C11a (hoxc11a)

| | |
|--------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-EP685369DIL-B |
| Abbreviation | hoxc11a |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q58EL0 |
| Product Type | Recombinant Protein |
| Immunogen Species | Danio rerio (Zebrafish) (Brachydanio rerio) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MFNSVNLGNF CSQTRKDRTS EFGDRTGCAS NIYLPSCITYY VPEFSAVSSF LPQGPSRQIT YPYSTNLSQV QPVRDVSYGL DPSSKWHHRS NYASCYSGED LVHRDCLPPS TMTEMLMKNE SVYSHHHQHH PNSSHPSTGF YTGVGKNNVL PQGFDRFFET AYCSTDNQSE HCLQKSEMNK LETSSQQPTA VSAAREPEKD PEDEEHTNS GSCTSAATKD GNASKSSHSS TPRTRKKRCP YSKFQIRELE REFFFNVIYIN KEKRLQLSRM LNLTDQRVVKI WFNQRRMKEK KLSRDRLQYF SGNPLL |
| Source | E.coli |
| Target Names | hoxc11a |
| Protein Names | Recommended name: Homeobox protein Hox-C11a |
| Expression Region | 1-306 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |