



Recombinant Human RUN domain-containing protein 3A (RUNDC3A)

Product Code	CSB-EP685437HU-B
Abbreviation	RUNDC3A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q59EK9
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MEASFVQTTM ALGLSSKKAS SRNVAVERKN LITVCRFSVK TLLEKYTAEP IDDSSEEFVN FAAILEQILS HRFKACAPAG PVSWFSSDGQ RGFWDYIRLA CSKVPNNCVS SIENMENIST ARAKGRWIR VALMEKRMSE YITTALRDTR TTRRFYDSGA IMLRDEATIL TGMLIGLSAI DFSFCLKGEV LDGKTPVVID YTPYLKFTQS YDYLTDEEER HSAESSTSED NSPEHPYLPL VTDEDSWYSK WHKMEQKFRI VYAQKGYLEE LVRLRESQLK DLEAENRRLQ LQLEEAAAQN QREKRELEGV ILELQEQLTG LIPSDHAPLA QGSKELTTPV LNQWPSLGTL NGAEGASNSK LYRRHSFMST EPLSAEASLS SDSQRLGEGT RDEEPWGPPIG KDPTPSMLGL CGSLASIPSC KSLASFKSNE CLVSDSPEGS PALSPS
Source	E.coli
Target Names	RUNDC3A
Protein Names	Recommended name: RUN domain-containing protein 3A Alternative name(s): Rap2-interacting protein 8 Short name= RPIP-8
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.