



Recombinant Danio rerio UPF0488 protein C8orf33 homolog (zgc:112185)

Product Code	CSB-BP685585DIL
Abbreviation	zgc:112185
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BJC4
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	MLEDEKITDA QSGEADPAIT PSSVASASAK KHKKKKKKKA GEGKDDQQRE TKTTEGETAK QETSELPDE QLSRELDWCI EQLELGLRTQ KTSSKQREEA SRALKTLRSS KAPLVKKRQV MRAISGDYRK KIEEEKNRQF KLIQSAMTSA RVTTVSECKP VFHRRERNT QPKNKTDGSQ ETQSSQETNE HTKTEDTKPF VFTKTHDFEC FDFDL
Source	Baculovirus
Target Names	zgc:112185
Protein Names	Recommended name: UPF0488 protein C8orf33 homolog
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.