



Recombinant Rat Keratin, type I cytoskeletal 18 (Krt18)

Product Code	CSB-BP685592RA
Abbreviation	Krt18
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BJY9
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	>85% (SDS-PAGE)
Sequence	SFTTRSTTF STNYRSLGSV RTPSQVRPA SSAASVYAGA GGSGSRISVS RSVWGGSVGS AGLAGMGGVQ TEKETMQDLN DRLASYLDKV KNLETENRRL ESKIREYLEK RGPQGV RDWG HYFKTIEDLR AQIFANSVDN ARIVLQIDNA RLAADDFRVK YETELAMRQS VESDIHGLRK VVDDTNITRL QLETEIEALK EELLFMKKNH EEEVQGLEAQ IASSGLTVEV DAPKSQDLSK IMADIRAQYE QLAQKNREEL DKYWSQQIEE STTVVTTKSA EIRDAETLL ELRRTLQTL EIDLDSMKNQN INLENNLGEV EARYRVQMEQ LNGVLLHLES ELAQTRAEGQ RQTQEYEA LL NIKVKLEAEI ATYRRILLEDG DDFSLNDALD SSNSMQTVQR TTTRKVVDGK VVSETNDTRV LRH
Source	Baculovirus
Target Names	Krt18
Protein Names	Recommended name: Keratin, type I cytoskeletal 18 Alternative name(s): Cytokeratin-18 Short name= CK-18 Keratin-18 Short name= K18
Expression Region	2-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	KRT18 encodes the type I intermediate filament chain keratin 18. Keratin 18, together with its filament partner keratin 8, are perhaps the most commonly found members of the intermediate filament gene family. They are expressed in single layer epithelial tissues of the body. Mutations in this gene have been linked to cryptogenic cirrhosis. Two transcript variants encoding the same protein have been found for this gene.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.