



Recombinant *Xenopus tropicalis* Vacuolar protein sorting-associated protein 26B (vps26b)

Product Code	CSB-EP685598XBF-B
Abbreviation	vps26b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5BKM4
Product Type	Recombinant Protein
Immunogen Species	<i>Xenopus tropicalis</i> (Western clawed frog) (<i>Silurana tropicalis</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSFFGFGPAA ELDIALTDGE SRRRVEHKTE DGKKEKYFLF YDGETVSGRV TVNLRNPGKR LEHQGLKIEF IGQIELYYDR GNHHEFVSLV KDLARPGEIS QSQSDFEFT HVEKPYESYT GQNVKLRVFL RATLSRRLND VVKEMDIVVH TLSTYPELNS SIKMEVGIED CLHIEFEYNK SKYHLKDVIV GKIYFLLVRI KIKHMEIDII KRETTGTGPN VYHENDTIAK YEIMDGAPVR GESIPIRLFL AGYELTPTMR DINKKFSVRY YLNLVLIDEE ERRYFKQQEV VLWRKGDIVR KSMHQAAIA SQRFEGTSHP ETRPQHSGAA AVEQEQE
Source	E.coli
Target Names	vps26b
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 26B Alternative name(s): Vesicle protein sorting 26B
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.