



# Recombinant Mouse Ras association domain-containing protein 5 (Rassf5)

<b>Product Code</b>	CSB-MP685744MO
<b>Abbreviation</b>	Rassf5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5EBH1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MASPAIGQRP YLLLLDPEPP RYLQSLGGTE PPPPARPRRC IPTALIPAAG ASEDRGGRRS GRRDPEPTPR DCRHARPVRP GLQPRLRLRP GSHRPRDVRS IFEQPQDPRV LAERGEGERHVELALRGGPG WCDLCGREVL RQALRCANCK FTCHSECRSL IQLDCRQKGG PALDRRSPES TLTPTLNQNV CKAVEETQHP PTIQEIKQKI DSYSNREKHC LGMKLSLEDGT YTGFIKVHLK LRRPVTVPAG IRPQSIYDAI KEVNPAATTD KRTSFYLPLD AIKQLHISST TTVSEVIQGL LKKFMVVDNP QKFALFKRIH KDGQVLFQKL SIADYPLYLR LLAGPDTDVL SFVLKENETG EVEWDAFSIP ELQNFLTILE KEEQDKIHQL QKKYNKFRQK LEEALRESQG KPG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	Rassf5
<b>Protein Names</b>	Recommended name: Ras association domain-containing protein 5 Alternative name(s): New ras effector 1 Regulator for cell adhesion and polarization enriched in lymphoid tissues Short name= RAPL
<b>Expression Region</b>	1-413
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This gene is a member of the Ras association domain family. It functions as a tumor suppressor, and is inactivated in a variety of cancers. The encoded protein localizes to centrosomes and microtubules, and associates with the GTP-activated forms of Ras, Rap1, and several other Ras-like small GTPases. The protein regulates lymphocyte adhesion and suppresses cell growth in response to activated Rap1 or Ras. Multiple transcript variants encoding different isoforms have been found for this gene.

**Reconstitution**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Shelf Life**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.