



Recombinant Chicken Stromal membrane-associated protein 2 (SMAP2)

Product Code	CSB-MP685765CH
Abbreviation	SMAP2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5F413
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MTGKSVRDVE RYQAVLGSL SEEENKYCAD CQAKGPRWAS WNIGVFICIR CAGIHRNLGV HISRVKSVNL DQWTQEIQIC MQEMGNGKAN RLYEAFLPEN FRRPQTDQAV EGFIRDKEYE KKYMDRSIDI NTFRKEKDDK WKKSNEPVSE RKLEPIVFEK VKMPQKKEET QQSRKSSPKS TEPVIDLLGL DAPVPSTLTN GRPCSLEKDL DLFASVGLNS DSRKVSGSMP TSGSAGSVPE NLNLFPEPGG KGEEAGKKQL SKDSILSLYG SQTPQLPAQG AMFMAPAQMA YPAAAYTSFP GVPPSSMMG GMAPSVGVM AQHGAAGMVT PMAIPAGYVG NVQAAVIGVP NGMMAAQQAG YVAGMAAVPQ PVYGVQPAQQ LQWNITQMTQ QMAGMNFYGA NGVMGYGQSM GGGGAQGSNQ SLSTQLWK
Source	Mammalian cell
Target Names	SMAP2
Protein Names	Recommended name: Stromal membrane-associated protein 2 Alternative name(s): Stromal membrane-associated protein 1-like
Expression Region	1-428
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.