



# Recombinant Chicken BAG family molecular chaperone regulator 5 (BAG5)

<b>Product Code</b>	CSB-YP685767CH
<b>Abbreviation</b>	BAG5
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5F486
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDMGNQHPSI KRLHEIQKEV KEIEQQVAVF SGLSTDRDYK KLERSLTKQL FEIDSVDTEG KGDIQQARKR AAQETERLLK ELEQNANHPR RLEIEAIFKE AQALVEREIT PFYQGGNCVN EEFEEGIQDV VLRLTQVKTG GKVSLRKARY RTLTKVCAVQ EIIESCAKRQ LSLPLSNDAAH PSVSKINSVM CEV NKARGTL IALLMGVSSN DTCRHLACVL TGLVADLDAL DVCGRTEIRN YRKEVVEEIN KLQKYLDLDE EANSTHAYDL AQNHSILKIE EIRKKLKEVN SLLLKTENAS DLYLGSKAEL QGLIAQLDEV SLGKNPCIRE ARRRAVIEVQ TLITYIDLKE ALGKRQMYAE QTAAEHQSHK AVWTVLGNLS QIQQEVISFD GNKTDKNYMR LEELLTKQLL ALDAVDPQGD ERCKAARKQA VKLAQNILYY LDMKTDEWEY
<b>Source</b>	Yeast
<b>Target Names</b>	BAG5
<b>Protein Names</b>	Recommended name: BAG family molecular chaperone regulator 5 Short name= BAG-5 Alternative name(s): Bcl-2-associated athanogene 5
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Target Details</b>	This protein is a member of the BAG1-related protein family. BAG1 is an anti-apoptotic protein that functions through interactions with a variety of cell apoptosis and growth related proteins including BCL-2, Raf-protein kinase, steroid hormone receptors, growth factor receptors and members of the heat shock protein 70 kDa family. This protein contains a BAG domain near the C-terminus, which could bind and inhibit the chaperone activity of Hsc70/Hsp70. Three transcript variants encoding two different isoforms have been found for this gene.



## Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.