



# Recombinant *Xenopus tropicalis* Target of rapamycin complex subunit lst8 (mlst8)

<b>Product Code</b>	CSB-YP686297XBF
<b>Abbreviation</b>	mlst8
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q5I0B4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Xenopus tropicalis</i> (Western clawed frog) ( <i>Silurana tropicalis</i> )
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNSNQGTVGS DPVILATAGY DHTVRFWQAH SGICTRTVQH QDSQVNSLEV TPDRSMIAAA GYQHIRMVYDL NSNNPNVIN YDGVSKNITS VGFHEDGRWM YTGGEDCMAR IWDLRSRNLQ CQRIFQVNAP INCVFLHPNQ AELIVGDQSG AIHIWDLKTD QNEQLIPETD VSINSVHIDP DASYMAAVNS SGNCFVWNLT GGLGEDLTQL IPKTKIPAHK RCALKCKFSP DSTLLATCSA DQTCKIWRTS NFSLMTELSI KSNPNGETSR GWMWDCAFSG DSQYIVTASS DNLARLWCVE TGEIKREYSG HQKAVVCLAF NDSVLG
<b>Source</b>	Yeast
<b>Target Names</b>	mlst8
<b>Protein Names</b>	Recommended name: Target of rapamycin complex subunit lst8 Short name= TORC subunit lst8 Alternative name(s): G protein beta subunit-like Short name= Gable Short name= Protein GbetaL MTOR associated protein, LST8 homolog
<b>Expression Region</b>	1-326
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.