



Recombinant Macaca thibetana Interleukin-6 (IL6)

Product Code	CSB-EP686313MDAO
Abbreviation	IL6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5I6E3
Product Type	Recombinant Protein
Immunogen Species	Macaca thibetana (Pere David's macaque) (Tibetan macaque)
Purity	>85% (SDS-PAGE)
Sequence	APV LPGEDSKDVA APHSQPLTSS ERIDKHIRYI LDGISALRKE TCNRSNMCES SKEALAENNL NLPKMAEKDG CFQSGFNEDT CLVKIITGLL EFEVYLEYLQ NRFESSEEQA RAVQMSTKAL IQFLQKKAKN LDAITTPPT TNASLLTKLQ AQNQWLQDMT THLILRSFKE FLQSSLRALR QM
Source	E.coli
Target Names	IL6
Protein Names	Recommended name: Interleukin-6 Short name= IL-6
Expression Region	28-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.