



Recombinant Human Krueppel-like factor 17 (KLF17)

Product Code	CSB-BP686385HU
Abbreviation	KLF17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q5JT82
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MYGRPQAEME QEAGELSRWQ AAHQAAQDNE NSAPILNMSS SSGSSGVHTS WNQGLPSIQH FPHSAEMLGS PLVSVEAPGQ NVNEGGPQFS MPLPERGMSY CPQATLTPSR MIYCQRMSPQ QQEMTIFSGP QLMPVGEPI PRVARPFGGN LRMPPNGLPV SASTGIPIMS HTGNPPVPYP GLSTVPSDET LLGPTVPSTE AQAFLPSMAQ MLPPQDAHDL GMPPAESQSL LVLGSQDSL SQPDSQEGPF LPEQPGPAPQ TVEKNSRPQE GTGRRGSSEA RPYCCNYENC GKAYTKRSHL VSHQRKHTGE RPYSCNWESC SWSFFRSDEL RRHMRVHTRY RPYKCDQCSR EFMRS DHLKQ HQBKTRPGPS DPQANNNGE QDSPPAAGP
Source	Baculovirus
Target Names	KLF17
Protein Names	Recommended name: Krueppel-like factor 17 Alternative name(s): Zinc finger protein 393
Expression Region	1-389
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.